

R-FORM

Fiberglass Cast Tape

R-FORM Features:

- Finished cast is strong and durable
- Radiolucent – superb to see through by X-ray
- Available in white only



Product Description	Product #
R-Form Casting Tape, Fiberglass - White 2" x 4yd (box=10 rolls)	4015-2, 10
R-Form Casting Tape, Fiberglass - White 3" x 4yd (box=10 rolls)	4015-3, 10
R-Form Casting Tape, Fiberglass - White 4" x 4yd (box=10 rolls)	4015-4, 10
R-Form Casting Tape, Fiberglass - White 5" x 4yd (box=10 rolls)	4015-5, 10

Casting Instructions:

1. Take the roll out of the package and immediately dunk into the water.
2. Make sure the water penetrates it thoroughly. Squeeze the roll hard while in the water to help penetrate and set off the resin within the water.
3. Take the roll out of the water and squeeze lightly the excess water. Key is to keep the roll wet.
4. USE A CUTTING STRIP.
5. When removing THE EZ CUTTER works best if the cast is less than 4-ply. Lay the blade to the side at a 45 degree angle.

- Tip #1:** The first layer of R-FORM Tape should be rolled gently. Then gently stretch second layer to release resin from casting tape.
- Tip #2:** When you feel the heat from the resin setting off, this is when you want to shape the cast to fit the patient's stump or limb.
- Tip #3:** Use HOT water for faster activation.
- Tip #4:** Use THE EZ CUTTER or a sharp pair of scissors when removing the cast.



ST&G Corporation USA
2691 Saturn Street
Brea, CA 92821

Telephone: (714)524-0663
FAX: (714)364-8113
www.stngco.com