



4-Hole Connector, Pyramid (With & Without Hole)

Instructions For Use

1217 & 1227T



1217A



1217S



1217T



1217T-B



1227T

1. Application & Description

- This adapter is exclusively for lower extremity prosthetic fittings.
- It connects any pyramid system to a 4-Hole adapter.
- The 1217 4-Hole Pyramid includes a through hole in the pyramid, while the 1227T does not.

2. Technical Specifications

P/N	Material	Build Height	Weight	Weight Limit
1217A	Aluminum	11mm	81g	158kg/350lbs
1217S	Stainless Steel	8mm	111g	136kg/300lbs
1217T	Titanium	8mm	70g	170kg/375lbs
1217T-B	Titanium	8mm	70g	170kg/375lbs
1227T	Titanium	8mm	73g	170kg/375lbs

- **Recommended Torque:** Always cross tighten Alignment Set Screws
- **M6 Flat Head Screw Torque:**
 - (Aluminum/ Titanium/ Stainless Steel): **10Nm**
- **Thread Locker:** Apply Loctite once alignment and connection are finalized.

3. Safety & Environment

- **Temperature Range:** Between -20°C to 80°C (-4°F to 176°F). 
- Rinse with clean water if exposed to salty, chlorinated, or otherwise potentially corrosive liquids.
- **Storage:** Store at room temperature; avoid wet, very hot, or extremely cold environments
- **Warranty:** 2 years, provided it is used as intended without modification, and only for the original intended user.
- Do not reuse on another prosthesis.
- **Conformity:** This product meets the requirements Council EU 2017/745 (MDR) for medical products. This product has been classified as a class I product according to the classification criteria outlined in appendix IX of the guidelines. Please keep this manual in safe place for future reference. 

1217_1227 IFU REV. B (7-21-2026)

Specifications subject to change.