



# 4-Hole Connector Receiver

1225 & 1215T

Instructions For Use



1225A



1225S/T



1225T-B



1215T

## 1. Application & Description

- This adapter is exclusively for lower extremity prosthetic fittings
- The 1225 is a Rotating Pyramid with hole
- The 1215T is Rotating Pyramid without hole
- **1225A is strictly for K2 ONLY.**
- Compatible with any 36mm threaded base, connecting a pyramid receiver system to a 36mm threaded receiver

## 2. Technical Specifications

| P/N     | Material            | Build Height | Weight | Weight Limit |
|---------|---------------------|--------------|--------|--------------|
| 1225A   | Aluminum            | 11mm         | 56g    | 158kg/350lbs |
| 1225S   | Stainless Steel     | 11mm         | 68g    | 136kg/300lbs |
| 1225T   | Titanium            | 11mm         | 39g    | 170kg/375lbs |
| 1225T-B | Titanium            | 11mm         | 39g    | 170kg/375lbs |
| 1215T   | Titanium w/out hole | 11mm         | 42g    | 170kg/375lbs |

- **Thread Locker:** Apply Loctite once alignment and connection are finalized

## 3. Safety & Environment

- **Temperature Range:** Between -20°C to 80°C (-4°F to 176°F).
- Rinse with clean water if exposed to salty, chlorinated, or otherwise potentially corrosive liquids.
- **Storage:** Store at room temperature; avoid wet, very hot, or extremely cold environments
- **Warranty:** 2 years, provided it is used as intended without modification, and only for the original intended user.
- Do not reuse on another prosthesis.
- **Conformity:** This product meets the requirements Council EU 2017/745 (MDR) for medical products. This product has been classified as a class I product according to the classification criteria outlined in appendix IX of the guidelines. Please keep this manual in safe place for future reference.



1225 1215T IFU REV. B (7-21-26)  
Specifications subject to change.