



Diagonal Offset 4-Hole Pyramid 1727

Instructions For Use



1727S/T



1727T-B

1. Application & Description

- This adapter is exclusively for lower extremity prosthetic fittings.
- It connects any 4-Hole to a pyramid receiver like a pylon pyramid receiver.
- The 1727 4-Hole Pyramid has a combination 5 mm lateral, and 7 mm diagonal offset to help correct misalignment.

2. Technical Specifications

P/N	Material	Build Height	Weight	Weight Limit
1727S	Stainless Steel	9mm	125g	136kg/300lbs
1727T	Titanium	9mm	92g	170kg/375lbs
1727T-B	Titanium	9mm	92g	170kg/375lbs

- **Recommended Torque:**
- **M6 Flat Head Screw Torque:**
 - (Titanium/ Stainless Steel): **10Nm**
- **Thread Locker:** Apply Loctite once alignment and connection are finalized.

3. Safety & Environment

- **Temperature Range:** Between -20°C to 80°C (-4°F to 176°F). 
- Rinse with clean water if exposed to salty, chlorinated, or otherwise potentially corrosive liquids.
- **Storage:** Store at room temperature; avoid wet, very hot, or extremely cold environments.
- **Warranty:** 2 years, provided it is used as intended without modification, and only for the original intended user.
- Do not reuse on another prosthesis. 
- **Conformity:** This product meets the requirements Council EU 2017/745 (MDR) for medical products. This product has been classified as a class I product according to the classification criteria outlined in appendix IX of the guidelines. Please keep this manual in safe place for future reference

1727 IFU REV. C (6-16-2026)

Specifications subject to change.