



Pylon, 30mm Aluminum Tube

Instructions For Use

1209A / 1210A

1209A



1210A



1. Application & Description

- This adapter is exclusively for lower extremity prosthetic fittings

2. Technical Specifications

P/N	Material	Pylon Length	Weight	Weight Limit
1209A	Aluminum	385mm	227g	170kg/375lbs
1210A	Aluminum	910mm	538g	136kg/300lbs

- **Compatibility:**

- The 385mm or 910mm Pylon, Al is designed for use with any 30mm tube clamp system.
- ST&G Tube clamps are highly recommended.

- **Installation Requirements:**

- **Shortening:** The tube must be shortened as necessary to fit the proper prosthesis height
- **Cutting:** End cut perpendicular to length of tube. Avoid using a tube cutter, as this will mushroom the ends causing potential binding within tube clamping surface
- **Finishing:** All burrs from both inner and outer surfaces must be removed properly
- **Fitment:** The tube must be completely inserted and in full contact with the inner base rim of the tube clamp. Never use a spacer in the tube clamp to make up length.
- **Fastening:** Tighten all screws that are side by side using a torque wrench

CAUTION: Do not clamp the tube in a vise or similar device, as this may damage the tube

3. Safety & Environment

- **Temperature Range:** Between -20°C to 80°C (-4°F to 176°F) 
- Rinse with clean water if exposed to salty, chlorinated, or otherwise potentially corrosive liquids
- **Storage:** Store at room temperature; avoid wet, very hot, or extremely cold environments
- **Warranty:** 2 years, provided it is used as intended without modification, and only for the original intended user
- Do not reuse on another prosthesis.
- **Conformity:** This product meets the requirements Council EU 2017/745 (MDR) for  cal products. This product has been classified as a class I product according to the classification criteria outlined in appendix IX of the guidelines. Please keep this manual in safe place for future reference.

1209A_1210A IFU REV. A (6-18-2026)

Specifications subject to change.