



# Double Ended Adapter

1201, 1202, 1203, & 1204

## Instructions For Use



1201T-B



1202T-B



1203T-B



1204T-B

### 1. Application & Description

- This adapter is for use in lower extremity prosthetic fittings.
- It connects two pyramid components using receiver connections on both ends.
- Four set screws secure the connection and allow for alignment adjustments.

### 2. Technical Specifications

| P/N     | Material | Build Height | Weight | Weight Limit |
|---------|----------|--------------|--------|--------------|
| 1201T-B | Titanium | 32mm         | 88g    | 170kg/375lbs |
| 1202T-B | Titanium | 45mm         | 99g    | 170kg/375lbs |
| 1203T-B | Titanium | 60mm         | 113g   | 170kg/375lbs |
| 1204T-B | Titanium | 75mm         | 126g   | 170kg/375lbs |

- **Recommended Torque:** Always cross tighten Alignment Set Screws
- **M8 Alignment Set Screw Torque:**
  - (Titanium/Stainless Steel): **15Nm**
- **Thread Locker:** Apply Loctite once alignment and connection are finalized

### 3. Safety & Environment

- **Temperature Range:** Between -20°C to 80°C (-4°F to 176°F).
- Rinse with clean water if exposed to salty, chlorinated, or otherwise potentially corrosive liquids.
- **Storage:** Store at room temperature; avoid wet, very hot, or extremely cold environments
- **Warranty:** 2 years, provided it is used as intended without modification, and only for the original intended user.
- Do not reuse on another prosthesis.
- **Conformity:** This product meets the requirements Council EU 2017/745 (MDR) for medical products. This product has been classified as a class I product according to the classification criteria outlined in appendix IX of the guidelines. Please keep this manual in safe place for future reference.



1201\_1202\_1203\_1204\_T-B IFU REV. A (10-5-26)

Specifications subject to change.