



Double Ended Adapter

1201, 1202, 1203, & 1204

Instructions For Use



1201T



1201S



1202T



1202S



1203T



1203S



1204T



1204S

1. Application & Description

- This adapter is for use in lower extremity prosthetic fittings.
- It connects two pyramid components using receiver connections on both ends.
- Four set screws secure the connection and allow for alignment adjustments.

2. Technical Specifications

P/N	Material	Build Height	Weight	Weight Limit
1201T	Titanium	32mm	88g	170kg/375lbs
1201S	Stainless Steel	32mm	150g	136kg/300lbs
1202T	Titanium	45mm	99g	170kg/375lbs
1202S	Stainless Steel	45mm	166g	136kg/300lbs
1203T	Titanium	60mm	113g	170kg/375lbs
1203S	Stainless Steel	60mm	192g	136kg/300lbs
1204T	Titanium	75mm	126g	170kg/375lbs
1204S	Stainless Steel	75mm	210g	136kg/300lbs

- **Recommended Torque:** Always cross tighten Alignment Set Screws
- **M8 Alignment Set Screw Torque:**
 - (Titanium/Stainless Steel): **15Nm**
- **Thread Locker:** Apply Loctite once alignment and connection are finalized

3. Safety & Environment

- **Temperature Range:** Between -20°C to 80°C (-4°F to 176°F).
- Rinse with clean water if exposed to salty, chlorinated, or otherwise potentially corrosive liquids.
- **Storage:** Store at room temperature; avoid wet, very hot, or extremely cold environments
- **Warranty:** 2 years, provided it is used as intended without modification, and only for the original intended user.
- Do not reuse on another prosthesis.
- **Conformity:** This product meets the requirements Council EU 2017/745 (MDR) for medical products. This product has been classified as a class I product according to the classification criteria outlined in appendix IX of the guidelines. Please keep this manual in safe place for future reference.



1201_1202_1203_1204 IFU REV. C (10-5-26)

Specifications subject to change.