



Heavy Duty 4-Ear Connector Pyramid

Instructions For Use

3230T



3230T

1. Application & Description

- This adapter is exclusively for lower extremity prosthetic fittings up to **375lb** max weight limit.
- This is for Trans-Tibial fittings.
- Be sure to layup and reinforce for the expected forces associated with higher and dynamic loads.

2. Technical Specifications

P/N	Material	Build Height	Weight	Weight Limit
3230T	Titanium	22mm	83g	170kg/375lbs

- **System Size:** 30mm

3. Safety & Environment

- **Temperature Range:** Between -20°C to 80°C (-4°F to 176°F).
- Rinse with clean water if exposed to salty, chlorinated, or otherwise potentially corrosive liquids.
- **Storage:** Store at room temperature; avoid wet, very hot, or extremely cold environments
- **Warranty:** 2 years, provided it is used as intended without modification, and only for the original intended user.
- Do not reuse on another prosthesis.
- **Conformity:** This product meets the requirements Council EU 2017/745 (MDR) for medical products. This product has been classified as a class I product according to the classification criteria outlined in appendix IX of the guidelines. Please keep this manual in safe place for future reference.



4. Reinforcement Recommendations

- The clinician is responsible for determining the adequate layup and reinforcement based on the specific patient's usage, activity level, and expected loading.

3230T IFU REV.E (3-17-2026)

Specifications subject to change.