



1. Application & Description

- This adapter is exclusively for lower extremity prosthetic fittings.
- It connects any pyramid system to a tube clamp.
- The four set screws anchor the mating system and allow for angular alignment.
- 1206T, 1207T, 1205T are compatible with the 1212T-SH short tube clamp for up to 375lb rating

2. Technical Specifications

P/N	Material	Build Height	Weight	Weight Limit
1206S	Stainless Steel	45mm	96g	125kg/275lbs
1206T	Titanium	45mm	61g	125kg/375lbs
1206T-B	Titanium	45mm	61g	125kg/375lbs
1207S	Stainless Steel	60mm	112g	125kg/275lbs
1207T	Titanium	60mm	71g	125kg/375lbs
1207T-B	Titanium	60mm	71g	125kg/375lbs
1205T	Titanium	100mm	128g	170kg/375lbs
1205T-B	Titanium	100mm	128g	170kg/375lbs

- **Recommended Torque:** Always cross tighten Alignment Set Screws
- **M8 Alignment Set Screw Torque:**
 - (Titanium/Stainless Steel): **15Nm**
- **Thread Locker:** Apply Loctite once alignment and connection are finalized

3. Safety & Environment

- **Temperature Range:** Between -20°C to 80°C (-4°F to 176°F).
- Rinse with clean water if exposed to salty, chlorinated, or otherwise potentially corrosive liquids.
- **Storage:** Store at room temperature; avoid wet, very hot, or extremely cold environments
- **Warranty:** 2 years, provided it is used as intended without modification, and only for the original intended user.
- Do not reuse on another prosthesis.
- **Conformity:** This product meets the requirements Council EU 2017/745 (MDR) for medical products. This product has been classified as a class I product according to the classification criteria outlined in appendix IX of the guidelines. Please keep this manual in safe place for future reference.



1205_1206_1207 IFU REV. B (6-29-26)

Specifications subject to change.