

3020 / 3070 PARASIL

Transtibial / Transfemoral Liner

Blister Reduction Comparison
Silicone vs. Parasil Silicone



Parasil Liner offers a superior durable design for your active patient and is available in both PIN LOCKING and CUSHION.

The silky smooth Parasil Liner is designed for high performance with minimal maintenance. Made with platinum cured medical grade mineral oil embedded silicone with a Parylene® coating gives the Parasil the unique ability to repel water and lower the co-efficiency of friction. This hydrophobic property is easy to clean and sanitize leaving it odor free.

The Parylene® surface reduces surface friction between the liner/skin and liner/socket which will help in the prevention of skin shear for fragile skin and blister reduction. The Parasil offers proprietary donning and doffing strips to aid the patient in easier application/removal. It is suitable for moderate to high activity patients who require cushioning over sensitive areas in the socket. The Parasil liner's opaque pigment gives a natural "skin tone" to all ethnic skin types.



Transtibial

P/N	Description - LOCKING	P/N	Description - CUSHION
3020L-18	Parasil Locking Liner, 18	3020C-18	Parasil Cushion Liner, 18
3020L-20	Parasil Locking Liner, 20	3020C-20	Parasil Cushion Liner, 20
3020L-22	Parasil Locking Liner, 22	3020C-22	Parasil Cushion Liner, 22
3020L-23.5	Parasil Locking Liner, 23.5	3020C-23.5	Parasil Cushion Liner, 23.5
3020L-25	Parasil Locking Liner, 25	3020C-25	Parasil Cushion Liner, 25
3020L-26.5	Parasil Locking Liner, 26.5	3020C-26.5	Parasil Cushion Liner, 26.5
3020L-28	Parasil Locking Liner, 28	3020C-28	Parasil Cushion Liner, 28
3020L-30	Parasil Locking Liner, 30	3020C-30	Parasil Cushion Liner, 30
3020L-32	Parasil Locking Liner, 32	3020C-32	Parasil Cushion Liner, 32
3020L-34	Parasil Locking Liner, 34	3020C-34	Parasil Cushion Liner, 34
3020L-36	Parasil Locking Liner, 36	3020C-36	Parasil Cushion Liner, 36

Transfemoral


P/N	Description - LOCKING	P/N	Description - CUSHION
3070L-30	Parasil Locking TF Liner, 30	3070C-30	Parasil Cushion TF Liner, 30
3070L-32	Parasil Locking TF Liner, 32	3070C-32	Parasil Cushion TF Liner, 32
3070L-34	Parasil Locking TF Liner, 34	3070C-34	Parasil Cushion TF Liner, 34
3070L-36	Parasil Locking TF Liner, 36	3070C-36	Parasil Cushion TF Liner, 36
3070L-38	Parasil Locking TF Liner, 38	3070C-38	Parasil Cushion TF Liner, 38
3070L-40	Parasil Locking TF Liner, 40	3070C-40	Parasil Cushion TF Liner, 40
3070L-45	Parasil Locking TF Liner, 45	3070C-45	Parasil Cushion TF Liner, 45




VISIT US AT STNGCO.COM
2691 Saturn Street | Brea, CA 92821
Tel: (714) 524-0663 | Fax: (714) 364-8113

CE
Authorized Representative
MDSS GmbH
Schiffgraben 41 | 30175 Hannover, Germany
Tel: +49 511 6262 8630 | Fax: +49 511 6262 8633


HOW TO MEASURE Transtibial Liner




1. Measure the circumference of the residual limb at 4cm from the distal end
2. Choose the liner size, which corresponds to the measurement or one smaller.
3. Never choose a liner size, which is larger than the measurement




ADAPT



PARASIL



CUSHION




MAGLINER


BK LINERS

WALL AVERAGES


HOW TO MEASURE Transfemoral Liner




1. Measure the circumference of the residual limb at 4cm from the distal end
2. Choose the liner size, which corresponds to the measurement or one smaller.
3. Never choose a liner size, which is larger than the measurement




ADAPT



PARASIL



CUSHION



MAGLINER

AK LINERS

WALL AVERAGES