




ADAPT

Patient Level: Low to Moderate
P/N 3012

- The ADAPT liner(s) offer an enhanced stabilizing distal matrix with flexible fabric that allows for a snug fit of various contours without restriction
- Conforms very nicely while providing comfort
- Suitable for moderate to high activity patients
- The enhanced distal matrix is part of the fabric covering of the liner and improves distal elongation control of the liner for both locking and cushion liner applications.
- Six-month warranty
- **ADAPT** and **AGILE** have the same thickness profile, thus are interchangeable if needed.



P/N	Description - LOCKING
3012L-18	ADAPT Locking Liner, 18
3012L-20	ADAPT Locking Liner, 20
3012L-22	ADAPT Locking Liner, 22
3012L-23.5	ADAPT Locking Liner, 23.5
3012L-25	ADAPT Locking Liner, 25
3012L-26.5	ADAPT Locking Liner, 26.5
3012L-28	ADAPT Locking Liner, 28
3012L-30	ADAPT Locking Liner, 30
3012L-32	ADAPT Locking Liner, 32
3012L-34	ADAPT Locking Liner, 34
3012L-36	ADAPT Locking Liner, 36



P/N	Description - CUSHION
3012C-18	ADAPT Cushion Liner, 18
3012C-20	ADAPT Cushion Liner, 20
3012C-22	ADAPT Cushion Liner, 22
3012C-23.5	ADAPT Cushion Liner, 23.5
3012C-25	ADAPT Cushion Liner, 25
3012C-26.5	ADAPT Cushion Liner, 26.5
3012C-28	ADAPT Cushion Liner, 28
3012C-30	ADAPT Cushion Liner, 30
3012C-32	ADAPT Cushion Liner, 32
3012C-34	ADAPT Cushion Liner, 34
3012C-36	ADAPT Cushion Liner, 36

AGILE

Patient Level: Active
P/N 3022


- The AGILE liner(s) offer a full matrix which enhances control for higher activity and secure fit.
- Helps in stabilizing the limb, yet provides comfort
- Suitable for moderate to high activity patients
- The enhanced distal matrix is part of the fabric covering of the liner and improves distal elongation control of the liner for both locking and cushion liner applications.
- A full matrix outer fabric covering that controls full length elongation and is available in both locking and cushion liner applications.
- Six-month warranty
- **ADAPT** and **AGILE** have the same thickness profile, thus are interchangeable if needed.



P/N	Description - LOCKING
3022L-18	AGILE Locking Liner, 18
3022L-20	AGILE Locking Liner, 20
3022L-22	AGILE Locking Liner, 22
3022L-23.5	AGILE Locking Liner, 23.5
3022L-25	AGILE Locking Liner, 25
3022L-26.5	AGILE Locking Liner, 26.5
3022L-28	AGILE Locking Liner, 28
3022L-30	AGILE Locking Liner, 30
3022L-32	AGILE Locking Liner, 32
3022L-34	AGILE Locking Liner, 34
3022L-36	AGILE Locking Liner, 36



P/N	Description - CUSHION
3022C-18	AGILE Cushion Liner, 18
3022C-20	AGILE Cushion Liner, 20
3022C-22	AGILE Cushion Liner, 22
3022C-23.5	AGILE Cushion Liner, 23.5
3022C-25	AGILE Cushion Liner, 25
3022C-26.5	AGILE Cushion Liner, 26.5
3022C-28	AGILE Cushion Liner, 28
3022C-30	AGILE Cushion Liner, 30
3022C-32	AGILE Cushion Liner, 32
3022C-34	AGILE Cushion Liner, 34
3022C-36	AGILE Cushion Liner, 36



ADAPT Locking
AGILE Locking

ADAPT Cushion
AGILE Cushion

Note: When placing an order, please use the below formula:

ADAPT TT Locking Liner (3012L)- indicate size

Example, if you were ordering a size 28, the part number would be: 3012L-28

AGILE TT Locking Liner (3022L)- indicate size


Example, if you were ordering a size 32, the part number would be: 3022L-32

HOW TO MEASURE Transtibial Liner


1. Measure the circumference of the residual limb at 4cm from the distal end
2. Choose the liner size, which corresponds to the measurement or one smaller.
3. Never choose a liner size, which is larger than the measurement




ST&G Liner Thickness Comparison




ADAPT LOCKING
AGILE LOCKING



PARASIL LOCKING



ADAPT/PARASIL
CUSHION
AGILE/PARASIL
CUSHION



PREMIER/HYDRO
MAGLINER